

DIOPTER TO RADIUS CONVERSION CHART

Diopter	Radius	Diopter	Radius	Diopter	Radius	Diopter	Radius	Diopter	Radius
23.00	14.67	34.67	9.73	40.12	8.41	45.67	7.39	51.12	6.60
24.00	14.06	34.75	9.71	40.25	8.39	45.75	7.38	51.25	6.59
25.00	13.50	34.87	9.68	40.37	8.36	45.87	7.36	51.37	6.57
26.00	12.98	35.00	9.64	40.50	8.33	46.00	7.34	51.50	6.55
27.00	12.50	35.12	9.61	40.67	8.30	46.12	7.32	51.67	6.53
28.00	12.05	35.25	9.57	40.75	8.28	46.25	7.30	51.75	6.52
29.00	11.64	35.37	9.54	40.87	8.26	46.37	7.28	51.87	6.51
30.00	11.25	35.50	9.51	41.00	8.23	46.50	7.26	52.00	6.49
30.12	11.21	35.67	9.46	41.12	8.21	46.67	7.23	52.12	6.48
30.25	11.16	35.75	9.44	41.25	8.18	46.75	7.22	52.25	6.46
30.37	11.11	35.87	9.41	41.37	8.16	46.87	7.20	52.37	6.44
30.50	11.07	36.00	9.38	41.50	8.13	47.00	7.18	52.50	6.43
30.67	11.00	36.12	9.34	41.67	8.10	47.12	7.16	52.67	6.41
30.75	10.98	36.25	9.31	41.75	8.08	47.25	7.14	52.75	6.40
30.87	10.93	36.37	9.28	41.87	8.06	47.37	7.12	52.87	6.38
31.00	10.89	36.50	9.25	42.00	8.04	47.50	7.11	53.00	6.37
31.12	10.85	36.67	9.20	42.12	8.01	47.67	7.08	53.12	6.35
31.25	10.80	36.75	9.18	42.25	7.99	47.75	7.07	53.25	6.34
31.37	10.76	36.87	9.15	42.37	7.97	47.87	7.05	53.37	6.32
31.50	10.71	37.00	9.12	42.50	7.94	48.00	7.03	53.50	6.31
31.67	10.66	37.12	9.09	42.67	7.91	48.12	7.01	53.67	6.29
31.75	10.63	37.25	9.06	42.75	7.89	48.25	6.99	53.75	6.28
31.87	10.59	37.37	9.03	42.87	7.87	48.37	6.98	53.87	6.27
32.00	10.55	37.50	9.00	43.00	7.85	48.50	6.96	54.00	6.25
32.12	10.51	37.67	8.96	43.12	7.83	48.67	6.93	54.12	6.24
32.25	10.47	37.75	8.94	43.25	7.80	48.75	6.92	54.25	6.22
32.37	10.43	37.87	8.91	43.37	7.78	48.87	6.91	54.37	6.21
32.50	10.38	38.00	8.88	43.50	7.76	49.00	6.89	54.50	6.19
32.67	10.33	38.12	8.85	43.67	7.73	49.12	6.87	54.67	6.17
32.75	10.31	38.25	8.82	43.75	7.71	49.25	6.85	54.75	6.16
32.87	10.27	38.37	8.80	43.87	7.69	49.37	6.84	54.87	6.15
33.00	10.23	38.50	8.77	44.00	7.67	49.50	6.82	55.00	6.14
33.12	10.19	38.67	8.73	44.12	7.65	49.67	6.79	55.25	6.11
33.25	10.15	38.75	8.71	44.25	7.63	49.75	6.78	55.50	6.08
33.37	10.11	38.87	8.68	44.37	7.61	49.87	6.77	55.75	6.05
33.50	10.07	39.00	8.65	44.50	7.58	50.00	6.75	56.00	6.03
33.67	10.02	39.12	8.63	44.67	7.56	50.12	6.73	56.25	6.00
33.75	10.00	39.25	8.60	44.75	7.54	50.25	6.72	56.50	5.97
33.87	9.96	39.37	8.57	44.87	7.52	50.37	6.70	56.75	5.95
34.00	9.93	39.50	8.54	45.00	7.50	50.50	6.68	57.00	5.92
34.12	9.89	39.67	8.51	45.12	7.48	50.67	6.66	57.25	5.90
34.25	9.85	39.75	8.49	45.25	7.46	50.75	6.65	57.50	5.87
34.37	9.82	39.87	8.47	45.37	7.44	50.87	6.63	57.75	5.84
34.50	9.78	40.00	8.44	45.50	7.42	51.00	6.62	58.00	5.82