



AccuCone^KTM Keratoconus Fitting System


Acculens[®] has been designing and manufacturing lenses for over 40 years. Over such time Acculens has conducted multiple studies at the University of Colorado covering many types of corneal distortion from Keratoconus, and Post Lasik, to various types of Corneal Trauma. Acculens provides a multitude of custom designs from corneal lenses to completely custom scleral lenses. Our experience in the design and research on these types of corneas is unsurpassed and our consultation staff is always available.

Accucone KTM is a lens indicated for early to advanced state keratoconus and can easily be adjusted to a custom fit. Accucone K provides excellent vision correction for keratoconus and is designed as a corneal lens.

DIAMETER	POWER	MATERIAL
7.90mm—10.2mm	±20.00	Any Material

The Accucone KTM was developed to cover a wide range of keratoconus patients.

Acculens provides the practitioner the ability to adjust the optical zone and AEL as desired for the best fit possible.



AEL can be ordered in Standard, Flat, and Steep.

Acculens is committed to providing exceptional quality, customer service and industry innovation for our customers, their patients and the ophthalmic community.

