

BACK VERTEX DISTANCE CONVERSION CHART

-	12 mm	+	-	12 mm	+	-	12 mm	+
-3.87	4.00	+4.25	-9.37	10.50	+12.00	-18.00	23.00	+31.75
-4.00	4.25	+4.50	-9.75	11.00	+12.75	-18.37	23.50	+32.62
-4.25	4.50	+4.75	-10.12	11.50	+13.37	-18.62	24.00	+33.62
-4.50	4.75	+5.0	-10.50	12.00	+14.00	-18.87	24.50	+34.75
-4.75	5.00	+5.25	-10.87	12.50	+14.75	-19.25	25.00	+35.75
-5.00	5.25	+5.62	-11.25	13.00	+15.50	-19.50	25.50	+36.75
-5.12	5.50	+5.87	-11.62	13.50	+16.12	-19.87	26.00	+37.75
-5.37	5.75	+6.12	-12.00	14.00	+16.75	-20.12	26.50	+38.87
-5.62	6.00	+6.50	-12.17	14.50	+17.50	-20.37	27.00	+40.00
-5.75	6.25	+6.75	-12.75	15.00	+18.25	-20.75	27.50	+41.00
-6.00	6.50	+7.00	-13.00	15.50	+19.00	-21.00	28.00	+42.25
-6.25	6.75	+7.37	-13.50	16.00	+19.75	-21.25	28.50	+43.50
-6.50	7.00	+7.62	-13.75	16.50	+20.50	-21.50	29.00	+44.50
-6.62	7.25	+8.00	-14.12	17.00	+21.50	-21.87	29.50	+45.66
-6.87	7.50	+8.25	-14.50	17.50	+22.25	-22.12	30.00	+47.00
-7.12	7.75	+8.50	-14.75	18.00	+23.00	-22.37	30.50	+48.12
-7.25	8.00	+8.87	-15.12	18.50	+23.75	-22.62	31.00	+49.50
-7.50	8.25	+9.12	-15.50	19.00	+24.75	-22.87	31.50	+50.75
-7.75	8.50	+9.50	-15.87	19.50	+25.50	-23.12	32.00	+52.12
-7.87	8.75	+9.75	-16.12	20.00	+26.37	-23.37	32.50	+53.50
-8.12	9.00	+10.12	-16.50	20.50	+27.11	-23.62	33.00	+54.62
-8.37	9.25	+10.37	-16.75	21.00	+28.12	-23.87	33.50	+56.12
-8.50	9.50	+10.75	-17.12	21.50	+29.00	-24.12	34.00	+57.50
-8.75	9.75	+11.00	-17.37	22.00	+29.87	-24.50	34.50	+59.50
-8.87	10.00	+11.37	-17.75	22.50	+30.87	-24.75	35.00	+60.62



WWW.ACCULENS.COM